

orlvision GmbH • Gewerbestr. 17 • D-35633 Lahnau

orlvision GmbH
Gewerbestr. 17
D-35633 Lahnau
Germany

Tel. +49 (0) 64 41-67 92 98-0
Fax. +49 (0) 64 41-67 92 98-99

info@orlvision.de
www.orlvision.de

Amtsgericht Wetzlar: HRB 5834
Geschäftsführer: Alexander Boll

UST-IDNr.: DE266282363

Lahnau, 11.09.2024

TO WHOM IT MAY CONCERN

Statement regarding equipment material compatibility approval

Dear Partner,

I, the undersigned, Alexander Boll of orlvision GmbH, hereby declare that the flexible video rhino-laryngoscopes manufactured by orlvision GmbH have been tested for material compatibility in collaboration with the UV-Smart D60.

These devices have successfully passed performance tests, including the evaluation of imaging components, in accordance with our internal procedures and protocols.

We would like to highlight that the colour of the handle may change during treatment with the UV-Smart D60, and the plastic surfaces may be affected depending on the number of treatment cycles.

Therefore, we recommend an annual inspection of the devices at orlvision if the UV-Smart D60 is used regularly.

Furthermore, we would like to clarify that orlvision devices can be reprocessed using the UV-Smart method, but it is the user's responsibility to validate the reprocessing method.

UV-Smart will list the devices tested with the UV-Smart D60 (see appendix) and will include a note regarding potential discolouration and alterations to the plastic surfaces.

orlvision GmbH

Alexander Boll
General manager

